



OLAES[®] HEMOSTATIC BANDAGE

Life-Assist
Helping Heroes Save Lives 

The next generation of the world's most versatile trauma bandage with built-in hemostatic gauze.



**TACTICAL MEDICAL
SOLUTIONS[®]**

PHONE 864.224.0081
TOLL FREE 1.888.822.6331
FAX 864.224.0064
E-MAIL info@tacmedsolutions.com
WEB tacmedsolutions.com

1250 Harris Bridge Road, Anderson, SC 29621, USA



Tactical Medical Solutions[®] is a global leading developer and manufacturer of pre-hospital medical gear for the tactical and civilian world.

We currently have over 2,000 products available and a presence in 70 countries around the world. Here at home, we service 7 of the top 10 police departments, including the NYPD and LAPD, as well as major fire and EMS agencies.

We're proud to say that when customers think of TacMed[™], they think uncompromising - uncompromising commitment to delivering top-quality products, tactical knowledge, customized training, and reliability in the field.

OLAES® HEMOSTATIC BANDAGE

Life-Assist
Helping Heroes Save Lives

| PRODUCT NUMBER | |
|--|--|
| OLAES® Hemostatic Bandage - 4" | 85-0020 |
| OLAES® Hemostatic Bandage - 4" Flat Packed | 85-0024 |
| DIMENSIONS | |
| 4" VERSION | 7.5"L x 5.5"W x 2"D 19cm L x 14cm W x 5cm D |
| 4" FLAT PACKED VERSION | 7.5"L x 5.5"W x 2"D 19cm L x 14cm W x 5cm D |
| WEIGHT | 2.7oz / 76.5g |
| US PATENTS | 7,652,190 & 8,163,973 B2 |



FEATURES

- HEMOSTATIC GAUZE**
Stocked with 3"x12' of ChitoGauze® PRO hemostatic gauze, ideal for a variety of applications and wound surfaces.
- TRUE PRESSURE CUP**
Generates focused pressure on the wound site to assist occluding damaged vessels. It can also serve as a rigid eye shield for ocular injuries per CoTCCC guidelines.
- WOUND PAD**
Covers the injury site to assist with sustained bleeding control and keep the wound clear of debris. Its non-adherent compartment keeps gauze tucked neatly inside for rapid access during application.
- VELCRO® CONTROL STRIPS**
Built into the elastic wrap to prevent accidental unrolling during application.
- ELASTIC WRAP**
Secures the gauze and wound pad to the injury site and assists in maintaining consistent applied pressure.
- SECURING CLIP**
Locks the bandage into place to free up the user's hands to treat other injuries or casualties.
- EASY-TEAR POINTS**
Intuitive, rapid access to the sterile bandage is achieved through multiple tear points indicated by red arrows at each package corner.

The **OLAES® Hemostatic Bandage** combines the globally recognized OLAES® Modular Bandage with battle-tested HemCon® ChitoGauze® PRO, creating the most comprehensive trauma bandage for multiple injury profiles.

One of few CoTCCC recommended hemostatic dressings, ChitoGauze® PRO is X-Ray detectable and the included gauze is located inside the bandage's wound pad for rapid access at the point of injury.

Having both a pressure bandage and hemostatic gauze in a single package streamlines your treatment process and saves space in your kit for additional medical items.

YOUR EQUIPMENT SHOULD NEVER BE A COMPROMISE.

